

TO STUDY THE INTERFERENCE OF LIGHT BY NEWTON'S RING APPARATUS

AIM

To determine the wave length of sodium light by using Newton's ring.

APPARATUS USED

Travelling microscope (for Newton's ring), Sodium vapour lamp, A thin plano-convex lens of large focal length, A convex lens of short focus, an optically plane glass plate.

THEORY AND FORMULA USED

When a monochromatic parallel beam of light is incident on a thin air film enclosed between a plano-convex lens and a plane glass plate (shown in Figure), the two reflected rays – one from top and another from bottom surface of thin air film (ray 2 and ray 3 shown in Figure) will interfere. Ray 1 and ray 4 shown in Figure will not take part in interference due to low coherency of the sodium light source. Around the point of contact, the points having equal thickness of air film lie on circles concentric with the point of contact. Hence due to interference between the two reflected rays a series of alternate dark and bright concentric rings will be produced. These are called Newton's rings. The thickness of air film at the contact point is infinitesimally small. But due to phase change of π by reflection from rarer to denser medium (air to glass), the central ring will be dark.

THEORY AND FORMULA USED

The wavelength of sodium light which is used for the Newton's ring setup is given by the formula:

$$\lambda = \frac{D_{n+m}^2 - D_n^2}{4mR}$$

where


λ = wavelength of sodium light

D_{n+m} = diameter of $(n+m)$ th ring


D_n = diameter of n th ring

R = Radius of curvature of the convex surface of the plano-convex lens.

SCHEMATIC OF THE NEWTON'S RING SETUP



Newton's ring arrangement



Interference between light waves reflected from different surfaces of plano-convex lens and plane glass plate.

PROCEDURE


1. Open the given link in a window

http://lo-au.vlabs.ac.in/laser-optics/Newtons_Rings_Wavelength_of_light/experiment.html

2. Select the medium (air), light source (Sodium) and click the knob light ON .
3. Order no. of the rings can be varied by moving the knob of microscope position.
4. According to the theory, the center of the interference fringes should be dark.
5. Move the microscope in horizontal direction to the left side of the fringes. Fix up the cross- wire tangentially to the $(n+m)^{\text{th}}$ bright ring and note down the readings of both the main scale and vernier scale. Then the microscope is moved in horizontal direction to the right side and should be fixed up tangentially to the successive bright fringes up to the 1st ring. Write down the readings of all successive rings. Then the microscope is moved to the right side of the fringes and should be fixed up tangentially to the 1st bright fringe. Write down the readings of both scales. In this way, all the readings of successive bright fringes should be noted down up to $(n+m)^{\text{th}}$ bright ring.

NEWTON'S RING SETUP (VIRTUAL)

Link : http://lo-au.vlabs.ac.in/laser-optics/Newtons_Rings_Wavelength_of_light/experiment.html



The image displays a virtual Newton's Ring experiment setup. On the left, a 3D illustration shows a microscope positioned over a glass plate with a lens. A circular inset at the top left shows the resulting Newton's rings pattern, which is a series of concentric dark and light rings. A scale bar at the bottom left of the inset shows measurements from 0 to 0.5 mm. On the right, a control panel titled 'VARIABLES' contains the following elements:

- Select Medium:** A dropdown menu currently set to 'Air'. A red arrow points to this menu with a blue box labeled 'Select'.
- Light Source:** A dropdown menu currently set to 'Sodium'. A red arrow points to this menu with a blue box labeled 'Select'.
- LIGHT ON:** A button to turn the light on. A red arrow points to this button with a blue box labeled 'Click'.
- Radius Of Lense :** Set to 50 cm, with a slider control below it.
- Microscope Position:** A slider control.
- Microscope Focus :** Set to 7.5 cm, with a slider control below it.
- RESET:** A button to reset the experiment.
- Show result:** A checkbox that is currently unchecked.
- RESULT:** A section header at the bottom of the control panel.

NEWTON'S RING SETUP (VIRTUAL)

Select the order no. of the ring by varying the microscope position

Note down the reading from the scale bar

Vary Microscope position

VARIABLES

Select Medium:
Air

Light Source
Sodium

LIGHT OFF

Radius Of Lense : 50 cm

Microscope Position

Microscope Focus : 7.5 cm

RESET

The image shows a virtual Newton's Ring setup. On the left, a 3D model of a microscope is shown with a yellow circular light source on its stage. A red arrow points from a text box to the microscope's objective lens. Below the microscope is a scale bar with markings from 0 to 0.5. A red arrow points from a text box to the scale bar. On the right, a control panel titled 'VARIABLES' contains several settings: 'Select Medium' (Air), 'Light Source' (Sodium), a 'LIGHT OFF' button, 'Radius Of Lense : 50 cm', 'Microscope Position' (a slider with a red dot), 'Microscope Focus : 7.5 cm', and a 'RESET' button. A red arrow points from a text box to the 'Microscope Position' slider.

OBSERVATIONS

- (A) Radius of the plano-convex lens, $R = \dots\dots\dots$ cm (note down from the reading)
- (B) Find the vernier constant of microscope :
- (C) Determination of Diameter of the Newton's rings

Table - I

Obs. No.	Ring No. n	Microscope Reading (in cm)						Diameter of the rings, $D = L_T \sim R_T$ (cm)
		L.H.S reading			R.H.S reading			
		MSR	VSR	TOTAL (L_T)	MSR	VSR	TOTAL (R_T)	
1.								
2.								
3.								
4.								
5.								

OBSERVATIONS

Table - II

Ring No. n	Diameter D_n	D_n^2 (cm ²)	$D_{n+m}^2 - D_n^2$ (cm ²)	Mean $D_{n+m}^2 - D_n^2$ cm ²
1.				
2.				
3.				
4.				
5.				

CALCULATIONS

$$\lambda = \frac{D_{n+m}^2 - D_n^2}{4mR}$$

RESULT

The wavelength of the sodium light is Å

PRECAUTIONS

- (1) The plano-convex lens should be of large radius of curvature.
- (2) The microscope is always moved in the same direction to avoid back lash error.
- (3) Here in virtual platform we are setting are crosswire on a bright ring but in physical experimental setup crosswire should be focused on a dark ring tangentially.

References

http://lo-au.vlabs.ac.in/laser-optics/Newtons_Rings_Wavelength_of_light/experiment.html